## Coherent manipulation of Andreev states in an atomic contact

<u>C. Janvier<sup>1</sup></u>, L. Tosi<sup>2</sup>, L. Bretheau<sup>3</sup>, Ç. Ö. Girit<sup>4</sup>, M. Stern<sup>1</sup>, P. Bertet<sup>1</sup>, P. Joyez<sup>1</sup>, D. Vion<sup>1</sup>, D. Esteve<sup>1</sup>, M. F. Goffman<sup>1</sup>, H.Pothier<sup>1</sup>, C. Urbina<sup>1\*</sup>.

<sup>1</sup> Quantronics Group, Service de Physique de l'État Condensé (CNRS, UMR 3680), IRAMIS, CEA-Saclay, 91191 Gif-sur-Yvette, France

<sup>2</sup>Present address: Centro Atómico Bariloche and Instituto Balseiro, Comisión Nacional de Energía Atómica (CNEA), 8400 Bariloche, Argentina

<sup>3</sup>Present address: Department of Physics, Massachusetts Institute of Technology, Cambridge, Massachusetts 02215, USA

<sup>4</sup>*Present address: Collège de France, 11 place Marcelin Berthelot, 75005 Paris, France* 

Adresse email de l'orateur : camille .janvier@cea.fr

Two superconducting electrodes coupled through a single conduction channel form the simplest Josephson weak-link. The Josephson effect can be fully described in terms of two quasiparticle states localized at the link and having energies within the superconducting gap. In this work we perform an exhaustive exploration of these so called Andreev bound states, including their coherent manipulation. The measurement is done by coupling a superconducting atomic contact to a microwave resonator. With this standard circuit QED architecture, we achieve single-shot measurements and observe quantum jumps.




Figure : Two tone spectrum of the Andreev transition. The amplitude of resonator signal is shown as a function of the excitation frequency and the phase difference on the weak-link. In dashed black is shown a fit for a transmission coefficient of 99.895%